
BRAEBON

BRAEBON Medical
1-100 Schneider Road
Kanata, ON, Canada, K2K
1Y2
01.613.831.6690

EC REP

CEpartner4U BV
Esdoornlaan 13
3951 DB Maarn
The Netherlands
www.cepartner4u.eu

Pneumo Flow Sensor

Intended Use: The PneuMo Flow Sensor is intended for use in the measurement of patient airflow at the proximal airway of patients using a CPAP device. The PneuMo Flow Sensor permits the simultaneous recording of both airflow and CPAP pressure when the device is connected to a pressure transducer and data acquisition system. The PneuMo Flow Sensor Kit is a limited-life device; two single-use hydrophobic (safety) filters must always be used to prevent both the spread of contaminants and moisture damage to the pressure transducer.

Connection

Blue tubing faces the patient; we recommend a short piece of flexible tubing between the PneuMo Flow Sensor and the patient connection. Ensure the tubing is upright to prevent kinking and moisture buildup.

Warning

The PneuMo Flow Sensor is a limited-life device. The PneuMo Flow Sensor is not sterile when supplied from the manufacturer. It is a clean device. The PneuMo Flow Sensor must be sterilized prior to use.

Hydrophobic filters are single-use devices; always use two new hydrophobic filters to prevent both the spread of contaminants and moisture damage to the pressure transducer.

Disinfection:

1. The PneuMo Flow Sensor is NOT sterile when supplied from the manufacturer. It is a clean device. The PneuMo Flow Sensor must be sterilized prior to use.
2. To sterilize the PneuMo Flow Sensor, disconnect the single-use hydrophobic filters. The PneuMo Flow Sensor can be sterilized using ethylene oxide gas or liquid sterilants (do not use temperatures in excess of 62° C or 145° F - steam autoclaving is not recommended). Ensure that the PneuMo Flow Sensor is completely dry before using.

EU Notice: any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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