
BRAEBON

			
<p>BRAEBON Medical 1-100 Schneider Road Kanata, ON, Canada, K2K 1Y2 01.613.831.6690 Toll-free (N.A.): 1.888.462.4841</p>	<p>EC REP CEpartner4U BV Esdorlaan 13 3951 DB Maarn The Netherlands www.cepartner4u.eu</p>		

en Thermistor and Thermocouple User Guide — All Models

Intended Use

Braebon thermistor and thermocouple airflow sensors are intended for use during sleep disorders studies as a qualitative measure of respiratory airflow for recording onto a data acquisition system. A thermistor airflow sensor requires an external power source, such as a battery; a thermocouple generates its own signal without an external power source.

Warning: Disinfect, rinse and dry airflow sensors prior to use. Sensors marked as single use must not be cleaned and reused. Airflow sensors are for diagnostic purposes only and are not intended as an apnea (apnoea) monitor.

Caution: US Federal law restricts this device to sale by or on the order of a physician. Airflow sensors are sensitive instruments and must be handled with care. Never immerse the battery enclosure at any time.

Note: Refer to www.braebon.com for more information. Product specifications are subject to change without notice.

How to Use

Refer to your sleep system user guide for manufacturer’s recommended settings. Red 1.5 mm touchproof connector is positive; black is negative.

General recording system settings:

- Low frequency filter: 0.1 Hz to 0.3 Hz (TC=1 sec. or longer)
- High frequency filter: 5 Hz to 15 Hz
- 50 / 60 Hz notch filter: optional

- Sensitivity / Gain: 7.5 µV/mm or 15 000 gain (increase or decrease as needed)
- Place sensor on patient such that sensor tips do not touch the skin and are exposed to airflow from both nose and / or mouth. Adjust plastic slide behind head or under chin. Use tape to securely affix sensor wire to both sides of face.

Specifications

Operating Temperature: 5°C (41° F)-32° C (90° F); REF **8514**: 10° C (50° F)-35° C (95° F)
Storage Temperature: 0°C (32° F)-60° C (140° F); REF **8514**: 0° C (32° F)-50 ° C (122° F)
Operating/Storage Humidity: 10–95% Non-condensing; REF **8514** operating 30-95%.

Disinfection

Caution: Never steam autoclave or gas sterilize the devices. Do not immerse the battery enclosure at any time.
Improper disinfection will void the warranty.

- Wash the sensor and attaching wire with soapy water and rinse thoroughly. Do not wet the black box which contains the battery.
- Allow the sensor to dry completely before reusing.
- You may disinfect with an EPA-registered hospital disinfectant.

Warranty

BRAEBON MEDICAL CORPORATION warrants to the first consumer that the airflow sensor (the “Sensor”), when shipped in its original container, will be free from defective workmanship, performance and materials and agrees that it will, at its option, either repair the defect or replace the defective Sensor or part thereof at no charge to the purchaser for parts or labor for a time period of seven months from the date of purchase. The cTherm Model 0514 is guaranteed for a period of three months. The warranty described herein shall be the sole and exclusive warranty granted by BRAEBON MEDICAL CORPORATION and shall be the sole and exclusive remedy available to the purchaser. **Improper disinfection will void the warranty.** Use of the Sensor constitutes total and complete acceptance of this warranty. Correction of defects, in the manner and for the time period described herein, shall constitute complete fulfillment of all liabilities and responsibilities of BRAEBON MEDICAL CORPORATION to the purchaser with respect to the Sensor and shall constitute full satisfaction of all claims, whether based on contract, negligence, strict liability or otherwise. In no event shall BRAEBON MEDICAL CORPORATION be liable, or in any way responsible, for any loss of revenues or damage, direct, incidental, or consequential, including property damage, loss of profit, or personal injury resulting from the use or misuse of, or the inability to use this product. Nor shall BRAEBON MEDICAL CORPORATION be liable, or in any way responsible, for any damages or defects in the Sensor which were caused by abuse, misuse, tampering, neglect, incorrect battery type, or repairs or attempted repairs performed by anyone other than an authorized service provider.

EU Notice: any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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