**ru**

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### cTherm Wegwerp

**REF** 0518-24 

**BRAEBON Medical Corporation**  
1-100 Schneider Road  
Kanata, ON, Canada, K2K 1Y2  
01.613.831.6690

**CEpartner4U BV**  
Esdoorniaan 13  
3951 DB Mearn  
The Netherlands  
www.cepartner4u.eu

© 2022 BRAEBON® Medical Corporation

D.0518.40600.1.1 2022.11.29 Release version 1.1

**en**


## Intended Use

The BRAEBON 0518 Disposable cTherm is a thin adhesive backed thermal sensor designed to easily attach to any cannula used in monitoring a patient's airflow. The 0518 Disposable cTherm is intended for use during sleep disorders studies as a qualitative measure of respiratory airflow for recording into a data acquisition system. Use the BRAEBON 0510BS Ultima Breath Sensor Interface Cable to connect the BRAEBON 0518 Disposable cTherm to standard polygraph equipment.

**Warning:** Sensors marked as single use (disposable) must NOT be cleaned and reused. Airflow sensors are for diagnostic purposes only and are not intended for use as an apnea monitor.

**Caution:** US Federal law restricts this device to sale by or on the order of a physician. Airflow sensors are sensitive instruments and must be handled with care. Never immerse the battery enclosure (of the 0510BS) at any time.

### Nasal thermistors


Plastic tab

Adhesive

Oral thermistor

## Setup

1. Remove the BRAEBON 0518 Disposable cTherm from its package.
2. Remove the paper backing from the BRAEBON 0518 Disposable cTherm to expose the adhesive.
3. Place the nasal prongs of the cannula through the holes on the sensor, ensuring the adhesive is down and the oral thermistor element is in front of the cannula, as illustrated. The blue will be facing toward the patient.
4. Gently squeeze to ensure the sensor sticks to the cannula. Wrap the two adhesive tabs around the cannula.

0518 Disposable cTherm on the PureFlow Cannula. The oral thermistor is placed toward the inside of the oral chamber.

5. Insert the plastic tab into the BRAEBON 0510BS thermistor connector. Fit the metal snap through the hole on the BRAEBON 0518 Disposable cTherm and snap closed.

Insert the plastic tab into the 0510BS thermistor connector.

EU Notice: any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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