



# COMPARISON BETWEEN LEVEL III PORTABLE MONITORING SYSTEM AND PSG ON 44 INLAB PATIENTS

## INTRODUCTION

AASM published guidelines in 2007 for the use of portable monitoring (PM) systems in the diagnosis of OSA. Concern remains about the validity of PM measurements. PM devices usually provide modules for automatic and manual scoring. Evaluation of current levels of performance require direct comparison between PM (automatic and manual scoring) and PSG data.

## METHOD

Prospective, random selection of 44 patients (Table 1) - suspected, on clinical history & examination, of having OSA - simultaneous PSG and PM recording in a sleep laboratory – PM = Braebon MediByte Jr. (MBJr), (Level III device - oronasal cannula, thorax RIP belt, finger pulse oximeter, body position sensor in MBJr) - AASM guidelines 2007 for sleep staging and respiratory events used - PM data analyzed using Braebon's automatic scoring module (MBA) and data scored manually (MBM) by experienced RPSGT. PSG data scored by the same tech - comparisons performed between MBA, MBM, and PSG. Sensitivity and specificity scores calculated for the presence and absence of OSA and for all degrees of OSA. MBJr data also recalculated with PSG TST.

## RESULTS

Table 1

Demographics & Clinical Profile					
Sex	M: 29, F: 15				
Age	Aver: 47yr (17-75yr)				
BMI	Aver: 33kg/m <sup>2</sup> (23-64)				
EDS	5	ESS	10	URS	2
FTG	4	Snore	84%	W/Aps	59%

(EDS – Excessive daytime sleepiness score (VAS 0-10); ESS – Epworth; URS – unrefreshing sleep (Likert 1-5); FTG – Fatigue score (VAS 0-10))

Table 2.2

Summary Chart of Raw Score Data contd.						
Average of all records (n=44)						
Events	MB Auto		MB Manual		PSG	
	Total	Index	Total	Index	Total	Index
Centrals	0	0	2	0	10	2
Obstructives	68	10	70	10	19	4
Mixed	0	0	0	0	4	1
Hypopneas	38	6	45	7	71	15
Apneas & Hypopneas	106	15	117	17	104	22
Desaturations	96	14	96	14		
Snoring	2669	375	2662	375	1014	185

Table 4

PSG	MB Auto				MB Manual			
	NPA	Mild	Moderate	Severe	NPA	Mild	Moderate	Severe
NPA	13	1	0	0	12	2	0	0
Mild	6	5	2	0	6	6	1	0
Moderate	0	3	2	0	0	0	5	0
Severe	0	0	2	10	0	0	1	11

Table 5

PSG	MB Auto		MB Manual	
	Sensitivity	Specificity	Sensitivity	Specificity
NPA	80% [1:5]	93% [1:14]	80% [1:5]	86% [1:7]
Mild	38% [2:3]	87% [1:8]	46% [1:2]	94% [1:17]
Moderate	40% [3:5]	87% [1:8]	100% [0:10]	81% [1:5]
Severe	83% [1:6]	100% [0:10]	92% [1:13]	100% [0:10]

(Assumption – PSG is the accurate score; then [ ] = misidentifications.)

Table 2.1

Summary Chart of Raw Score Data						
Sum of all records (n=44)						
Events	MB Auto		MB Manual		PSG	
	Total	Index	Total	Index	Total	Index
Centrals	0	0	99	12	425	91
Obstructives (Aps)	2990	431	3068	448	847	188
Mixed	0	0	3	0	191	35
Hypopneas (Hyps)	1656	243	1993	292	3109	673
Aps & Hyps	4646	674	5146	752	4581	986
Desaturations	4208	610	4245	616		
Snoring	117425	16498	117107	16513	44612	8136

Table 3

(n = 44)	AHI	NPA (< 5)	(%)	Mild (>=5, <15)	(%)	Moderate (>=15, <30)	(%)	Severe (>=30)	(%)
MBJr Auto	15	19	43%	9	20%	6	14%	10	23%
MBJr Manual	17	18	41%	8	18%	7	16%	11	25%
PSG	22	14	32%	13	30%	5	11%	12	27%

(NPA = non-pathological apnea)

Table 6

	MB Auto	MB Manual	PSG
Total recording/sleep time	6.90	6.84	4.65
AHI by own time	15 [-32%]	17 [-23%]	22
AHI by PSG time	23 [+5%]	25 [+14%]	22

(Assumption – respiratory events probably occurred during sleep only, therefore, dividing by true sleep time will give more accurate AHI. [ ] = variance from PSG result.)

## CONCLUSIONS

- The automatic scoring system in this device is within the range of the sensitivities of PSG studies (75-88%).
- Having the data reviewed manually significantly improves the sensitivities of the OSA severity classes.
- Identifying TST improves the concordance substantially.
- Clinical decisions that seem warranted from this study based on MBJr alone, that is, using automatic scoring module data only, include:
  - If the MBJr identifies OSA, good agreement with PSG [13:14]; meets minimum AHI criteria for diagnosis; treatment justified.
  - If MBJr indicates NPA, cannot rely on it [wrong 1:5 patients], non-treatment unjustified.
  - If MBJr indicates mild or moderate OSA and therefore, possible oral appliance treatment, not reliable [wrong 1:6 patients], OA treatment not justified; if OA therapy not a consideration, CPAP treatment justified.
  - If MBJr indicates severe OSA, CPAP treatment justified [100% agreement].

AUTHOR AND INSTITUTION

Support: Nil

Dr. John B. Carlile  
MB, BCh, DPM, PhD, FRCPS  
drcarlile@gmail.com

Medical Director  
Limestone City Sleep Laboratory  
Kingston, Ontario, Canada

Acknowledgements:  
Thanks to night techs who performed the studies  
and to Kris Marks who scored all the records.



LIMESTONE CITY  
SLEEP LABORATORY