<b>BRAEBON</b>					
<span></span>					
<b>Canada:</b> <div>1-100 Schneider Road</div> <div>Kanata, ON, Canada, K2K 1Y2</div> <div>01.613.831.6690</div> <div>Toll-free (N.A.): 1.888.462.4841</div>	<div><div><span><span></span></span></div><div><span><span></span></span></div></div>	<div><div><span><span></span></span></div><div><span><span></span></span></div></div>	<div><div><span><span></span></span></div><div><span><span></span></span></div></div>	<div><div><span><span></span></span></div><div><span><span></span></span></div></div>	<div><div><span><span></span></span></div><div><span><span></span></span></div></div>
	<b>EC</b>	<b>REP</b>	<b>CEpartner4U BV</b>	<b>CE</b>	
			<b>Esdoomrtaan 13</b>		
			<b>3951 DB Maarn</b>		
			<b>The Netherlands</b>		
			<b>www.cepartner4u.eu</b>		

**EN** **Ultima Pressure Transducer User Guide – All Models**

**Intended Use**

Ultima Pressure Transducers are for use during sleep disorder studies as a measure of respiratory airflow (hypopneas and apneas) and/ or esophageal pressure for recording onto a data acquisition system. The battery-powered sensor uses either a disposable cannula with a 0,2-micron hydrophobic filter, or an air-filled disposable esophageal catheter. Attach either cannula or catheter to the patient, connect to the Pressure Transducer input, and plug directly into the electrode headbox or a DC amplifier.

**Warning:** Pressure Transducers are for diagnostic purposes only and are not intended as an apnea (apnoea) monitor. Always use either a new disposable air-filled Esophageal Catheter (BRAEBON part #0591) with the necessary safety filter (part #0583Pes) or a single-patient use BRAEBON Airflow Canula with the necessary safety filter (part #0583) to prevent the spread of contaminants between patients and to prevent damage to the pressure sensor. To prevent dust contamination to the pressure sensor always keep safety filters attached to the unit and replace immediately prior to the next patient. Use a brand new cannula and filter after each patient test.

**Caution:** US Federal law restricts this device to sale by or on the order of a physician. Ensure correct battery type is used or damage to the Ultima Airflow Pressure Transducer will result. Do not mix battery types and do not insert the batteries backwards or damage may occur to the pressure transducer. Do not immerse the pressure transducer in any liquid. Do not steam autoclave or gas sterilize the pressure sensor or damage will result.

**Note:** Product specifications are subject to change without notice. On Dual Pressure Transducer use esophageal catheter on oral input port only.

**How to Use**

Positive input ports (+) may be used for either cannula or esophageal catheter connection. Typically, the oral input (+) port on the Dual Pressure Transducer is used for esophageal catheter recording. Negative input port (-) may be used to connect to a CPAP circuit.

**Note:** Refer to <http://www.braebon.com/faqs.html> for more information about input / output connections and the Model 0580Pes esophageal pressure transducer calibration.

Refer to your sleep system user guide for manufacturer’s recommended amplifier settings. Red 1.5 mm touchproof output connector is positive; black is negative.

General recording system settings:

- Low frequency filter: 0.1 Hz (TC=1 sec or longer)
- High frequency filter: 5 Hz to 15 Hz
- 50 / 60 Hz notch filter: optional

- Sensitivity / Gain: 2 mV/cm or 4,000 gain (increase or decrease as needed)

To subtract a known pressure from the nasal baseline signal, for example during CPAP titration, connect the CPAP circuit to the BRAEBON Safety Filter (part #0583) and connect the filter to the nasal differential port on the Ultima Dual Pressure Sensor. Note that you will only be able to affect the baseline on the nasal output signals.

The Ultima Dual Pressure Sensor contains four 20-turn potentiometers which allow for precise adjustment of both baseline and gain settings for the nasal and oral outputs. The baseline is set to zero volts in the factory and should rarely need adjusting. To increase gain, turn the gain screw clockwise; to decrease gain, turn the gain counterlockwise.

Two AA alkaline batteries will last approximately 30 to 50 8-hour recordings (i.e., 30-50 nights). The pressure transducer interface will automatically turn off after twenty hours of continuous use.

**Specifications**

Operating Temperature: 5°C (41°F)-32°C (90°F)

Storage Temperature: 0°C (32°F)-60°C (140°F)

Operating/Storage Humidity: 10–95% Non-condensing

**Disinfection and Cleaning**

**Warning:** BRAEBON airflow cannulas are single-patient use only and are not to be used by more than one patient.

The Ultima Pressure Transducer does not contact the patient; wipe with an isopropyl alcohol pad as needed. For additional information, refer to the APIC guidelines for selection and use of disinfectants (American Journal of Infection Control, Vol. 18, number 2, April, 1990, p. 99-117). **Steam autoclaving is not recommended.**

**Warranty**

BRAEBON MEDICAL CORPORATION warrants to the first consumer that the Ultima Pressure Transducer (the "Sensor"), when shipped in its original container, will be free from defective workmanship, performance and materials and agrees that it will, at its option, either repair the defect or replace the defective Sensor or part thereof at no charge to the purchaser for parts or labor for a time period of one year from the date of purchase. The warranty described herein shall be the sole and exclusive warranty granted by BRAEBON MEDICAL CORPORATION and shall be the sole and exclusive remedy available to the purchaser. Use of the Sensor constitutes total and complete acceptance of this warranty. Correction of defects, in the manner and for the time period described herein, shall constitute complete fulfillment of all liabilities and responsibilities of BRAEBON MEDICAL CORPORATION to the purchaser with respect to the Sensor and shall constitute full satisfaction of all claims, whether based on contract, negligence, strict liability or otherwise. In no event shall BRAEBON MEDICAL CORPORATION be liable, or in any way responsible, for any loss of revenues or damage, direct, incidental, or consequential, including property damage, loss of profit, or personal injury resulting from the use or misuse of, or the inability to use this product. Nor shall BRAEBON MEDICAL CORPORATION be liable, or in any way responsible, for any damages or defects in the Sensor which were caused by abuse, misuse, tampering, neglect, incorrect battery type, or repairs

or attempted repairs performed by anyone other than an authorized service.

**EU Notice:** any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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