

MediByte™ Snoring Report

Pre-treatment Summary Report (baseline – no treatment)


Patient Name:	Smith, Mary		Patient ID:	7732767863
Date of Birth:	03/14/1954		Chart Code:	0987654321
Weight:	135.0 lbs		Study Date:	07/12/07
Height:	5' 2"	157 cm	Age:	55
BMI:	24.8		Sex:	Female
Waist: 34"	Hip: 32"	Waist-Hip Ratio: 1.06	Referring Physician:	Dr. Snooze

Comments: An oral appliance with follow-up MediByte test is recommended.

Total Recording Time(TRT) : 365.9 minutes

Respiratory and Snoring Events	Total#	Index	Duration (sec.)			Cardiac					
			Mean	Min.	Max.	Avg HR:	70.4	Min HR:	57.0	Max HR:	98.0
Central Apneas	0	0.0	0.0	0.0	0.0	Oximetry					
						Mean SpO2		93.9%			
Obstructive Apneas	70	11.5	36.1	10.2	118.1	Min SpO2		88.0%			
						Max SpO2		98.0%			
Mixed Apneas	0	0.0	0.0	0.0	0.0	SpO2 Range		%	Minutes		
Hypopneas	14	2.3	22.4	11.0	56.6	90-100 %		99.8%	364.8		
Apnea+Hypopnea	84	13.8	33.8	10.2	118.1	80-89 %		0.2%	0.6		
Snoring	2005	328.8	0.8	0.2	2.4	70-79 %		0.0%	0.0		
Desaturations	40	6.6	31.9	8.3	80.0	60-69%		0.0%	0.0		
RERAs	0	0.0	0.0	0.0	0.0	50-59%		0.0%	0.0		
						< 50 %		0.0%	0.0		

Snoring Volume Table dB(A)						Body Position	Supine	Prone	Left Side	Right Side	Total Non Supine
Body Position	Supine	Non-Supine	Right Side	Left Side	Prone	% Time in Pos	73.9%	0.0%	12.4%	13.6%	26.1%
Min	50.0	50.0	50.0	50.0	-	Snoring events	1679	0	263	63	326
Max	99.9	95.4	85.4	95.4	-	Apnea+ Hypopnea	79	0	2	3	5
Mean	68.5	60.9	58.2	63.4	-	Apnea- Hypopnea Index.	17.5	0.0	2.6	3.6	3.1



Patient	Apnea-Hypopnea Index	Severe	Moderate	Mild	Normal
Mary Smith	13.8	>30	15 to 30	5 to 15	<5

*Respiratory events are defined in the Assisted Scoring User Settings and in the User Guide. Final clinical decisions and degree of accuracy are the sole responsibility of the clinician using this software. Powered by BRAEBON®